STERILE EO

Fiber Optic Light Probe

DIRECTIONS FOR USE

(English)

Circulating Nurse:

1. Peel open the outer pouch letting the inner pouch drop on to the sterile field use caution not to touch the inner pouch.

Scrub Nurse:

2. Peel open the sterile pouch, place illuminated pick onto the sterile field turning tray up-side-down and dislodging the pick handle by pressing on the tray with thumb.
3. Break and remove paper tape from fiber optic coil.
4. Hand off the adaptor end of the cable to the circulating nurse.

Circulating Nurse:

5. Insert the adaptor end of the cable into the appropriate light source (using adaptors as needed).*

**** Please see back cover of DFU for product numbers of different style light probes and for appropriate adaptors.***

WARNING

1. This sterile disposable medical device should not be reused. After use, all parts of this disposable pack should be properly disposed of.
2. If any component in this pack is received in a defective condition, do not use and notify Escalon immediately.

Please Direct Inquiries to:

Escalon Medical Corp.
2440 South 179th Street
New Berlin, WI 53146 USA
Tel 800-433-8197
262-821-9182
Fax 262-821-9927

Authorized European Representative

Drew Scientific Limited Co.
Sowerby Woods Industrial Estate
Barrow in Furness
Cumbria LA14 4QR
United Kingdom
Tel: +44 1229 432089
Fax: +44 1229 432096

LIGHT PROBE PRODUCTS

Product Number		Description	Description of Use
REF CE ₀₄₇₃	TR9800 -11 -61	Standard Light Probe, 20 ga	Intraocular Illumination
REF CE ₀₄₇₃	TR9801 -11 -61	45° Light Pick, 20 ga	Intraocular Illumination with Tissue Manipulation
REF CE ₀₄₇₃	TR9806 -11 -61	Retinal Knife, 20 ga	Intraocular Illumination with Tissue Manipulation
REF	TR9809 -11	Hardten Iris Retractor, 23 ga	Hooked tip for manipulation and retraction of Iris for optimal visualization of anterior vitreous.
REF CE ₀₄₇₃	TR9821 -11 -61	Wide Angle Light Probe	Wide View Intraocular Illumination
REF	TR9879 -11	Siepsner Shielded Light Probe	Identifying Debris and Fibers in the Corneal Interface or epithelium.
REF	TR9880 -11	DCR Illuminator	Transillumination of Tissue

U.S. Patent No. 4,878,487

-11 Series Probes	Alcon® Clinitex® CooperVision® Escalon® Heslin Mackool® Miva® OMS® Storz® Trek®	-61 Series Probes	Escalon® Trek® Grieshaber®
Connect (with appropriate adaptor) to the following light sources:		Connect directly to the following light source:	

ADAPTORS FOR -11 SERIES

Please contact Escalon Customer Service for assistance in determining proper adaptor.

REF	TR9811-2	Silver	Alcon®, Heslin Mackool®, Mira®, Storz® and others
REF	TR9814	Blue	OMS® and Clinitex®
REF	TR9815	Red	Escalon® and Trek®

Please Direct Inquiries to:

Escalon Medical Corp.
2440 South 179th Street
New Berlin, WI 53146 USA
Tel 800-433-8197
262-821-9182
Fax 262-821-9927

Authorized European Representative

Drew Scientific Limited Co.
Sowerby Woods Industrial Estate
Barrow in Furness
Cumbria LA14 4QR
United Kingdom
Tel: +44 1229 432089
Fax: +44 1229 432096

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